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IVD

REF:

EIA-1468

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01/05

Legal Manufacturer:

DRG Instruments GmbH, Germany
Division of DRG International, Inc
Frauenbergstr. 18, D-35039 Marburg
Telefon: +49 (0)6421-17000
Fax: +49-(0)6421-1700 50
Internet: www.drg-diagnostics.de
E-mail: drg@drg-diagnostics.de

Distributed by:

DRG International, Inc.
1167 U.S. Highway 22 E
Mountainside, NJ 07092 USA
Telephone: (908) 233-2079
Fax: (908) 233-0758
Internet: www.drg-international.com
E-mail: corp@drg-international.com

CE 0197

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Legal Manufacturer:

DRG Instruments GmbH, Germany
Division of DRG International, Inc
Frauenbergstr. 18, D-35039 Marburg
Telefon: +49 (0)6421-17000
Fax: +49-(0)6421-1700 50
Internet: www.drg-diagnostics.de
E-mail: drg@drg-diagnostics.de

Distributed by:

DRG International, Inc.
1167 U.S. Highway 22 E
Mountainside, NJ 07092 USA
Telephone: (908) 233-2079
Fax: (908) 233-0758
Internet: www.drg-international.com
E-mail: corp@drg-international.com

CE 0197

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| | | |
|----|------------------------------------|---|
| 1 | 简介(INTRODUCTION)..... | 2 |
| 2 | 测试原理(PRINCIPLE OF THE TEST)..... | 2 |
| 3 | 注意事项(PRECAUTIONS)..... | 2 |
| 4 | 试剂盒组成(KIT COMPONENTS)..... | 3 |
| 5 | 样品(SPECIMEN)..... | 4 |
| 6 | 实验步骤(TEST PROCEDURE)..... | 4 |
| 7 | 期望值(EXPECTED VALUES)..... | 6 |
| 8 | 质量控制 (Quality Control) | 7 |
| 9 | 实验技术指标(ASSAY CHARACTERISTICS)..... | 7 |
| 10 | 使用注意事项(LIMITATIONS OF USE)..... | 8 |
| 11 | 相关法律事项(LLEGAL ASPECTS)..... | 8 |
| 12 | 参考文献(REFERENCES)..... | 9 |

Symbols used with DRG ELISA's

| Symbol | English | Deutsch | 中文 | Français | Espanol | Italiano | Portugues | Dansk | Svenska | Ελληνικά |
|---|---------------------|------------------------------|---------|-----------------------------------|-----------------------------|------------------------------|--------------------------------------|----------------------------|----------------------------|-----------------------------|
|  | European Conformity | CE-Konfirmitätskennzeichnung | CE 认证 | Conformité aux normes européennes | Conformidad europea | Conformità europea | Conformidade com as normas europeias | Europæisk overensstemmelse | Europeisk överensstämmelse | Ευρωπαϊκή Συμμόρφωση |
|  | In vitro diagnostic | In-vitro-Diagnostikum | 体外诊断试剂 | Diagnostic in vitro | Diagnóstico in vitro | Diagnostica in vitro | Diagnóstico in vitro | In vitro diagnostik | Diagnostik in vitro | in vitro διαγνωστικό |
|  | Catalogue number | Katalog-Nr. | 产品目录号 | Référence | No de catálogo | No. di Cat. | Catálogo n.º | Katalognummer | Katalognummer | Αριθμός καταλόγου |
|  | Lot. No. | Chargen-Nr. | 产品批号 | No. de lot | Número de lote | Lotto no | No do lote | Lot nummer | Batch-nummer | Αριθμός Παρτίδος |
|  | Storage Temperature | Lagerungstemperatur | 贮藏温度 | Temperature de conservation | Temperatura de conservación | Temperatura di conservazione | Temperatura de conservação | Opbevaringstemperatur | Förvaringstemperatur | Θερμοκρασία αποθήκευσης |
|  | Expiration Date | Mindesthaltbarkeits-datum | 失效日期 | Date limite d'utilisation | Fecha de caducidad | Data di scadenza | Prazo de validade | Udløbsdato | Bäst före datum | Ημερομηνία λήξης |
|  | Legal Manufacturer | Hersteller | 法定生产厂商 | Fabricant | Fabricante | Fabbricante | Fabricante | Producent | Tillverkare | Κατασκευαστής |
| Distributed by | Distributor | Distributeur | 代理经销商 | Distributeur | Distribuidor | Distributore | | | | |
|  | User's Manual | Arbeitsanleitung | 使用手册 | Mode d'emploi | Instrucciones de empleo | Istruzioni d'uso | Instruções de uso | Brugermanual | Användar manual | Εγχειρίδιο χρήστη |
| Content | Content | Inhalt | 内容 (组成) | Conditionnement | Contenido | Contenuto | Conteúdo | Indhold | Innehåll | Περιεχόμενο |
| Volume/No. | Volume / No. | Volumen/Anzahl | 数量 | Volume/Quantité | Volumen/Número | Volume/Quantità | Volume/Número | Volumen/antal | Volym/antal | Όγκος/αριθ.. |
| Microtiter wells | Microtiterwells | Mikrotiterwells | 微孔 | Barrettes de microtitration | Pocillos de la Microplaca | Pozzetti della micropiastra | Alvéolos de microtitulação | Mikrotiterbrønne | Brunnar i Mikrotiterplatta | Πηγαδάκια Μικροπιλοδοτήσεως |
| Antiserum | Antiserum | Antiserum | 抗血清 | Antisérum | Antisero | Antisiero | Anti-soro | Antiserum | Antiserum | Αντιπώρος |
| Enzyme Conjugate | Enzyme Conjugate | Enzym Konjugat | 酶联物 | Conjugué enzymatique | Conjugado enzimático | Tracciante enzimatico | Conjugado enzimático | Enzymkonjugat | Enzymkonjugat | Συζευγμένο ενζύμο |
| Enzyme Complex | Enzyme Complex | Enzym Komplex | 酶复合物 | Complex enzymatique | Complex enzimático | Complesso enzimatico | Complexo enzimático | Enzymkomplex | Enzymkomplex | Σύμπλοκο ενζύμου |
| Substrate Solution | Substrate Solution | Substrat Lösung | 底物液 | Solution substrat | Solución de sustrato | Soluzione di substrato | Solução de substrato | Substratopløsning | Substratlösning | Διάλυμα υποστρώματος |
| Stop Solution | Stop Solution | Stopplösung | 终止液 | Solution d'arret | Solución de paro | Soluzione di arresto | Solução de paragem | Stopopløsning | Stopp lösning | Διάλυμα τερματισμού |
| Zero Standard | Zero Standard | Nullstandard | 0 标准液 | Standard 0 | Standard 0 | Standard 0 | Padrão zero | Standard 0 | Standard 0 | Πρότυπο Μηδέν |

| Standard | Standard | Standard | 标准液 | Standard | Calibrador | Calibratore | Calibrador | Standard | Standard | Πρότυπα |
|-------------------|-------------------|---------------------------|---------|--------------------|--------------------|---------------------------------|--------------------|--------------|---------------|-----------------------------|
| Control | Control | Kontrolle | 对照 (质控) | Controle | Control | Controllo | Controlo | Kontrol | Kontroll | Έλεγχος |
| Assay Buffer | Assay Buffer | Assay Puffer | 实验缓冲液 | Tampon d'essai | Tampón de ensayo | Tampone del test | Tampão de teste | Assay buffer | Assay Buffer | Ρυθμιστικό Διάλυμα Εξέτασης |
| Wash Solution | Wash Solution | Waschlösung | 清洗液 | Solution de lavage | Solución de lavado | Soluzione di lavaggio | Solução de lavagem | Vaskebuffer | Tvätt lösning | Διάλυμα πλύσεως |
| 1N NaOH | 1N NaOH | 1N NaOH | 1N NaOH | 1N NaOH | 1N NaOH | 1N NaOH (idrossido di sodio 1N) | 1N NaOH | 1N NaOH | 1N NaOH | 1N NaOH |
| Sample Diluent | Sample Diluent | Probenverdünnungsmedium | 样品稀释液 | | | | | | | |
| Conjugate Diluent | Conjugate Diluent | Konjugatverdünnungsmedium | 酶联物稀释液 | | | | | | | |

1. 简介 (Introduction) :

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-
- 7.** Abelev, G.I., Alpha-fetoprotein in ontogenesis and its association with malignant tumors. *Adv. Cancer Res.* 14: 295, 1971.
 - 8.** Scardino, P. T., Cox, H. D., Waldemann, T. A., McIntire, K. R., Mittermeyer, B. and Javadpour, N., The value of the serum tumor markers in the staging and prognosis of germ cell tumors of the testis. *J. Urol.* 118:994, 1977.
 - 9.** Bosl, G.J., Lange, P. H., Fraley, E. E., Goldman, A., Nochomovitz, L. E., Rosai, J., et al., Human chorionic gonadotropin and alpha-fetoprotein in the staging of nonseminomatous testicular cancer. *Cancer* 47:328, 1981
 - 10.** Lange, P.H., McIntire, K.R. and Waldmann, T. A., Hakala, T.R. and Fraley, E. E., Serum alpha-fetoprotein and human chorionic gonadotropin in the diagnosis and management of nonseminomatous germ-cell testicular cancer. *Medical Intelligence* 259:1237, 1976.
 - 11.** Javadpour, N., McIntire, K. R. and Waldmann, T. A., Human chorionic gonadotropin (HCG) and alpha-fetoprotein (AFP) in sera and tumor cells of patients with testicular seminoma. *Cancer* 42:2768, 1978.
 - 12.** Waldmann, T.A. and McIntire, K. R., Serum Alpha-Fetoprotein Levels in Patients with Ataxia-Telangiectasia. *Lancet* 2:1112, 1972.
 - 13.** Belanger, L., Tyrosinémie Héritaire et Alpha-1-Fetoprotéine. II. Recherche Tissulaire Comparée de l'Alpha-Foetoprotéine dans Deux Cas de Tyrosinémie Héritaire. Considerations sur l'Ontogenèse de la Foetoprotéine Humain (French). *Pathol. Bio.*21:457, 1973.
 - 14.** Kew, M. C., Purves, L.R. and Bersohn, I., Serum Alpha-Fetoprotein Levels in Acute Viral Hepatitis. *Gut* 14:939, 1973.
 - 15.** Endo, Y., Kanai, K., Oda, T., Mitamura, K., Iino, S. and Suzuki, H., Clinical Significance of alpha-Fetoprotein in Hepatitis and Liver Cirrhosis. *Ann. N. Y. Acad. Sci.* 259:234, 1975.
 - 16.** Purves, L. R. and Purves M., Serum Alpha-Fetoprotein. VI. The Radioimmunoassay Evidence for the Presence of AFP in the Serum of Normal People and During Pregnancy. *S. Afr. Med. J.* 46:1290, 1972.
 - 17.** Sepalla, M. and Ruoslahti, E., Alpha-Fetoprotein: Physiology and Pathology During Pregnancy and Application to Antenatal Diagnosis. *J. Perinat. Med.* 1:104, 1973.
 - 18.** Engvall, E., *Methods in Enzymology*, Volume 70, Van Vunakis, H. and Langone, J.J. (eds.) Academic Press, New York, NY, 419 (1980).
 - 19.** Uotila, M., Ruoslahti, E. and Engvall, E., *J. Immunol. Methods*, 42, 11 (1981).

ii.