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CA 125

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IVD

REF:

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EIA-3939

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01/05

Legal Manufacturer:

DRG Instruments GmbH, Germany
Division of DRG International, Inc
Frauenbergstr. 18, D-35039 Marburg
Telefon: +49 (0)6421-17000
Fax: +49-(0)6421-1700 50
Internet: www.drg-diagnostics.de
E-mail: drg@drg-diagnostics.de

Distributed by:



DRG International, Inc.
1167 U.S. Highway 22 E
Mountainside, NJ 07092 USA
Telephone: (908) 233-2079
Fax: (908) 233-0758
Internet: www.drg-international.com
E-mail: corp@drg-international.com

CE 0197

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1	简介(INTRODUCTION).....	2
2	测试原理(PRINCIPLE OF THE TEST).....	2
3	注意事项(PRECAUTIONS).....	2
4	试剂盒组成(KIT COMPONENTS).....	3
5	样品(SPECIMEN).....	4
6	实验步骤(TEST PROCEDURE).....	5
7	期望值(EXPECTED VALUES).....	6
8	实验技术指标(ASSAY CHARACTERISTICS).....	6
9	使用注意事项(LIMITATIONS OF USE).....	8
10	相关法律事项(LLEGAL ASPECTS).....	8
11	参考文献(REFERENCES).....	9

Symbols used with DRG ELISA's

Symbol	English	Deutsch	中文	Français	Espanol	Italiano	Portugues	Dansk	Svenska	Ελληνικά
	European Conformity	CE-Konformitätskennzeichnung	CE 认证	Conformité aux normes européennes	Conformidad europea	Conformità europea	Conformidade com as normas europeias	Europaeisk overensstemmelse	Europeisk överensstämmelse	Ευρωπαϊκή Συμμόρφωση
	In vitro diagnostic	In-vitro-Diagnostikum	体外诊断试剂	Diagnostic in vitro	Diagnóstico in vitro	Diagnostica in vitro	Diagnóstico in vitro	In vitro diagnostik	Diagnostik in vitro	in vitro διαγνωστικό
	Catalogue number	Katalog-Nr.	产品目录号	Référence	No de catálogo	No. di Cat.	Catálogo n.º	Katalognummer	Katalognummer	Αριθμός καταλόγου
	Lot. No.	Chargen-Nr.	产品批号	No. de lot	Número de lote	Lotto no	No do lote	Lot nummer	Batch-nummer	Αριθμός Παρτίδος
	Storage Temperature	Lagerungstemperatur	贮藏温度	Temperature de conservation	Temperatura de conservación	Temperatura di conservazione	Temperatura de conservação	Opbevaringstemperatur	Förvaringstemperatur	Θερμοκρασία αποθήκευσης
	Expiration Date	Mindesthaltbarkeits-datum	失效日期	Date limite d'utilisation	Fecha de caducidad	Data di scadenza	Prazo de validade	Udløbsdato	Bäst före datum	Ημερομηνία λήξης
	Legal Manufacturer	Hersteller	法定生产厂商	Fabricant	Fabricante	Fabbricante	Fabricante	Producent	Tillverkare	Κατασκευαστής
Distributed by	Distributor	Distributeur	代理经销商	Distributeur	Distribuidor	Distributore				
	User's Manual	Arbeitsanleitung	使用手册	Mode d'emploi	Instrucciones de empleo	Istruzioni d'uso	Instruções de uso	Brugermanual	Användarmanual	Εγχειρίδιο χρήστη
Content	Content	Inhalt	内容 (组成)	Conditionnement	Contenido	Contenuto	Conteúdo	Indhold	Innehåll	Περιεχόμενο
Volume/No.	Volume / No.	Volumen/Anzahl	数量	Volume/Quantité	Volumen/Número	Volume/Quantità	Volume/Número	Volumen/antal	Volym/antal	Όγκος/αριθ..
Microtiter wells	Microtiterwells	Mikrotiterwells	微孔	Barrettes de microtitration	Pocillos de la Microplaca	Pozzetti della micropiastra	Alvéolos de microtitulação	Mikrotiterbrønne	Brunnar i Mikrotiterplatta	Πηγαδάκια Μικροπιλοδοτήσεως
Antiserum	Antiserum	Antiserum	抗血清	Antisérum	Antisero	Antisiero	Anti-soro	Antiserum	Antiserum	Αντιορός
Enzyme Conjugate	Enzyme Conjugate	Enzym Konjugat	酶联物	Conjugué enzymatique	Conjugado enzimático	Tracciante enzimatico	Conjugado enzimático	Enzymkonjugat	Enzymkonjugat	Συζευγμένο ενζύμο
Enzyme Complex	Enzyme Complex	Enzym Komplex	酶复合物	Complex enzymatique	Complex enzimático	Complesso enzimatico	Complexo enzimático	Enzymkomplex	Enzymkomplex	Σύμπλοκο ενζύμου
Substrate Solution	Substrate Solution	Substrat Lösung	底物液	Solution substrat	Solución de sustrato	Soluzione di substrato	Solução de substrato	Substratopløsning	Substratlösning	Διάλυμα υποστρώματος
Stop Solution	Stop Solution	Stopplösung	终止液	Solution d'arret	Solución de paro	Soluzione di arresto	Solução de paragem	Stopopløsning	Stopp lösning	Διάλυμα τερματισμού
Zero Standard	Zero Standard	Nullstandard	0 标准液	Standard 0	Standard 0	Standard 0	Padrão zero	Standard 0	Standard 0	Πρότυπο Μηδέν

DRG CA 125 ELISA EIA-3939

Standard	Standard	Standard	标准液	Standard	Calibrador	Calibratore	Calibrador	Standard	Standard	Πρότυπα
Control	Control	Kontrolle	对照 (质控)	Controle	Control	Controllo	Controlo	Kontrol	Kontroll	Έλεγχος
Assay Buffer	Assay Buffer	Assay Puffer	实验缓冲液	Tampon d'essai	Tampón de ensayo	Tampone del test	Tampão de teste	Assay buffer	Assay Buffer	Ρυθμιστικό Διάλυμα Εξέτασης
Wash Solution	Wash Solution	Waschlösung	清洗液	Solution de lavage	Solución de lavado	Soluzione di lavaggio	Solução de lavagem	Vaskebuffer	Tvätt lösning	Διάλυμα πλύσεως
1N NaOH	1N NaOH	1N NaOH	1N NaOH	1N NaOH	1N NaOH	1N NaOH (idrossido di sodio 1N)	1N NaOH	1N NaOH	1N NaOH	1N NaOH
Sample Diluent	Sample Diluent	Probenverdünnungsmedium	样品稀释液							
Conjugate Diluent	Conjugate Diluent	Konjugatverdünnungsmedium	酶联物稀释液							

1. 简介 (Introduction) :

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8. Davis, H.M., Zurawski, V.R. Jr., Bast, R.C. Jr., and Klug, T.L.: Characterization of the CA 125 antigen associated with human epithelial ovarian carcinomas. *Cancer Res.* 46:6143-6148, 1986.
9. Berek, J.S., Knapp, R.C., Malkasian, G.D., Lavin, P.T., Whitney, C., Niloff, J.M. and Bast, R.C. Jr.: CA 125 serum levels correlated with second-look operations among ovarian cancer patients. *Obstet. Gynecol.* 67:685-689, 1986.
10. Niloff, J.M., Bast, R. C., Jr., Schaetzel, E.M., and Knapp, R.C.: Predictive value of CA 125 antigen levels at second-look procedures in ovarian cancer. *Am. J. Obstet. Gynecol.* 151:981-986, 1985.
11. Atack, D.B., Nisker, J.A., Allen, H.H., Tustanoff, E.R. and Levin, L.: CA 125 surveillance and second-look laparotomy in ovarian carcinoma. *Am. J. Obstet. Gynecol.* 154:287-289, 1986.
12. Niloff, J.M., Klug, T.L., Schaetzel, E., Zurawski, V.R., Jr., Knapp, R.C., and Bast, R.C., Jr.: Elevation of serum CA 125 in carcinomas of the fallopian tube, endometrium, and endocervix. *Am. J. Ob. Gynecol.* 148:1057-1058, 1984.
13. Ruibal, A., Encabo, G., Martinez-Miralles, E., Murci, C., Capdenba B.A. Salgado, and Martinex-Vasquez, J.M.: CA 125 serum levels in non-malignant pathologies. In *Bull. Cancer (Paris)*, Masson Press 71:145-148, 1984.
14. Eerdekens, M.W., Nouwen, E.J., Pollet, D.E., Briers, T.W. and DeBroe, M.E.: Placental alkaline phosphatase and cancer antigen 125 in sera of patients with benign and malignant diseases. *Clin. Chem.* 31:687-690, 1985.
15. Halila, H., Stenman, U.H. and Seppala, M.: Ovarian cancer antigen CA 125 in pelvic inflammatory disease and pregnancy. *Cancer* 57:1327-1329, 1986.
16. Barbieri, R.L., Niloff, J.M., Bast, R.C., Jr., Schaetzel, E. Kistner, R.W. and Knapp, R.C.: Elevated serum concentrations of CA 125 in patients with advanced endometriosis. *Fert. and Steril.* 45:630-634, 1986.
17. Ricolleau, G., Chatal, J.F, Fumoleau, P., Kremer, M., Douillard, J. Y, and Curtet, C.: Radiomunoassay of the CA 125 antigen in ovarian carcinomas: advantages compared with CA 19-9 and CEA. *Tumor Biol.* 5:151-159, 1984.
18. Canney, P.A., Moore, M., Wilkinson, P.M., and James, R.D.: Ovarian cancer antigen CA 125: a prospective clinical assessment of its role as a tumor marker. *Br. J. Cancer* 50:765-769, 1984.
19. Krebs, H.B., Goplerud, D.R., Kilpatrick, S.J., Myers, M.B. and Hunt, A.: Role of CA 125 as tumor marker in ovarian carcinoma. *Obstet. Gynecol.* 67:473-477, 1986.
20. Kivinen, S., Kuoppala, T., Leppilampi, M., Vuori, J. and Kauppila, A: Tumor-associated antigen CA 125 before and during the treatment of ovarian carcinoma. *Obstet. Gynecol.* 67:468-472, 1986.
21. Einhorn, N., Bast, R.C., Jr., Knapp, R.C., Tjernberg, B. and Zurawski, V.R., Jr.: Preoperative evaluation of serum CA 125 levels in patients with primary epithelial ovarian cancer. *Obstet. & Gynecol.* 67:414-416, 1986.
22. Zurawski, V.R., Broderick, S.F., Pickens, P., Knapp, R.C., and Bast, R.C. Jr., Serum CA 125 levels in a Group of Non-hospitalized Women: Relevance for the early detection of Ovarian Cancer. *Obstet. & Gynecol.* 69:606-611, 1987.
23. Niloff, J.M., Knapp, R.C., Schaetzel, E.Reynolds, C., and Bast, R.C., Jr., Jr.: CA 125 levels in obstetric and gynecologic patients. *Obstet. & Gynecol.* 64:703-707, 1984.
24. Bast, R. C., Jr., Klug, T.L., Schaetzel, E., Lavin, O., Niloff, J.M., Greber, T.F., Zurawski, V.R., Jr., and Knapp, R.C.: Monitoring human ovarian carcinoma with a combination of CA 125, CA 19 9, and CEA. *Am. J. Obstet. Gynecol.* 149:553-559,1984.
25. Muto, M.G., Finkler, N.J., Kassis, A.I., Lepisto, E.M., and Knapp, R.C.: Human antimurine antibody responses in ovarian cancer patients and undergoing radioimmunotherapy with the murine monoclonal antibody OC 125. *Gynecol. Oncol.* 38:244-248, 1990.
26. Turpeinen, U., lehtovirta, P., Alfthan, H., and Stenman, U.: Interface by human antiidiotypic antibodies not neutralized by mouse IgG but removed by chromatography. *Clin. Chem.* 36(7):1333-1338, 1990.
27. Reinsberg, J., Heydweiller, A., Wagner, U., Pfell, K., Oehr, P., Krebs, D.: Evidence for interaction of human anti-idiotype antibodies with the CA 125 determination in a patient after radiomunodetection. *Clin. Chem.* 36(1):164-167, 1990.
28. Reinsberg, J. and Gast, B.: Human antiidiotypic antibodies against OC 125 strongly interfere with one-step assays of CA 125 employing OC 125 and M II antibodies. *Clin. Chem.* 40(6):951-952, 1994.

29. Kenemans, P., van Kamp, G.J., Oehr, P. and Verstraeten, R.A.: Heterologous double-determinant immunoradiometric assay Ca 125 II: reliable second-generation immunoassay for determining CA 125 in serum. Clin. Chem. 39:2509-13, 1993.